

OUR NEWEST PRODUCT DEVELOPMENT
THE BEST OF THE BESTEST™

NEW!!!!!!!



GrabHold™ Stone Baskets Nitinol and Tipless design

PERFORMANCE

MED-Fibers products including the **GrabHold™** Stone Baskets matches the highest standards on all existing product specifications available on the world wide market. We are the **BEST OF THE BESTEST™**

RELIABILITY

MED Fibers Laser Fibers and accessories like the **GrabHold™** Stone Basket are GUARANTEED compatible for safe and effective use with all FDA approved ureteroscopes. Our products are registered under:

FDA 510(k): K 124003

We also have ISO 9001:2008

and ISO 13485: 2003

and also CE 0459

QUALITY

MED-Fibers GrabHold™ Stone Baskets and Laser Fiber delivery systems are manufactured to provide optimum performance and 100% quality inspected before leaving production.

QUALITY ENFORCEMENT IS OUR MOTTO, BECAUSE WE ARE THE BEST OF THE BESTEST™.



MED-Fibers, Inc. with its key management and the combined optical glass, and optical fiber experience is your partner for the next generation of technology including accessories as a part of our product portfolio. The permanent performing product improvements will provide you with the quality you expect from your vendor and we will be the **BEST** partner for **CUSTOMER SERVICE RESPONSE TIME, LEAD TIME AND PRODUCT QUALITY.** Our technicians are selected and trained on an ongoing basis for the understanding of work performance, product handling and the importance of **Customer Satisfaction.** We as a team are on a 24/7 stand by for all of your requests and as your support team to provide you with solutions. **WE ARE:**

MED-Fiber, Inc. GrabHold™ Stone Baskets with its Nitinol and tipless design allows to position the basket directly against the mucosal lining to reduce bleeding, maintain a great and clear working field and reduce trauma.

The memory **characteristics** of the nitinol wire allows the basket to return to its originally shape when the extractor is developed.

The **flexibility** of the baskets sheets allows you to fully deflect the ureteroscope during stone manipulation.

The **scissor style** handle is designed to keep the basket and stone within the focal range of the ureteroscope.

Available products:

- 1.7FR 4WIRE 120CM NITINOL- TIPLESS
- 1.7FR 4WIRE 120CM NITINOL- TIPLESS, flexible
- 1.9FR 4WIRE 120CM NITINOL- TIPLESS
- 1.9FR 4WIRE 120CM NITINOL- TIPLESS, flexible
- 2.2FR 4WIRE 120CM NITINOL- TIPLESS
- 2.2FR 4WIRE 120CM NITINOL- TIPLESS, flexible
- 2.4FR 4WIRE 120CM NITINOL- TIPLESS,
- 2.4FR 4WIRE 120CM NITINOL- TIPLESS, flexible

ALL THE PRODUCTS ARE SINGLE USE AND STERILE

MED-Fibers, Inc. provides you with the highest standard of laser fiber devices and accessories like our **GrabHold™** Stone Baskets with affordable prices and the shortest delivery times on the market. This includes custom made and OEM labeled products.

QUALITY ENFORCEMENT IS OUR MOTTO AND WE GUARANTEE YOU RECEIVE THE BEST AVAILABLE PRODUCT WITH THE MOST COMPETITIVE PRICES.



THE BEST OF THE BESTEST™





DESIGN & ENGINEERING

MED-Fibers **GrabHold™** Stone Baskets and Laser Delivery Devices are designed and engineered by highly educated and trained engineers and support staff.

MANUFACTURING


MED-Fibers **GrabHold™** Stone Baskets and Laser Delivery Devices are produced in State of the Art facilities , with 0% tolerance of failure rates. We are the **BEST OF THE BESTEST™**

TURNKEY SOLUTIONS

MED-Fibers experienced Team with huge and advanced optical fiber and surgical accessory knowledge will support you in all of your projects and assists you within 48 hours response time. Your and our idea and our quick product solution, can help you to introduce your next product quicker into the market. For more information on any of our products or services please visit us on the Web at:

www.MED-Fibers.com

Service Features and Benefits



MED Fibers, Inc Laser fiber delivery systems and accessories are 100% tested for internal validation of the assembly performance. This gives us the opportunity to verify our 0% tolerance failure rate. We are the **BEST OF THE BESTEST™**

For questions, please contact our technical department or sales support team which will provide you with satisfied answers.



SERVICES AVAILABLE

Technical Support

Phone: +1 (480) 598 4090 ext. 311

Customer Support / Clinical Support

Phone: +1 (480) 598 4090 ext. 301

General Support

Phone: +1 (480) 598 4090 ext. 303

Quality Department

Phone: +1 (480) 598 4090 ext. 304

Financial Department

Phone: +1 (480) 598 4090 ext. 302



THE BEST OF THE BESTEST™



Headquarter:

7404 W. Detroit St, # 140
Chandler, AZ 85226, USA
Phone: +1 (480) 598 4090
Fax: +1 (480) 940 1049
Email: info@med-fibers.com